neodisher® FA

Alkaline detergent for the automated reprocessing of thermostable and thermolabile instruments
- liquid concentrate -

Main fields of application:

Automated cleaning of surgical instruments, anaesthetic equipment, theatre shoes and other medical utensils as well as baby bottles in hospitals and medical practices.

Automated cleaning of laboratory glassware in medical laboratories.

Characteristics:

neodisher FA is a versatile cleaning agent which is very gentle on material having a good cleaning power in respect of blood, protein, culture media and gelatine. neodisher FA is free of oxidation agents and surfactants.

neodisher FA is suitable for surgical instruments, anaesthetic equipment, laboratory glassware or plastics. Anodised aluminium surfaces must be tested for resistance. With titanium and titanium alloys, there might be colour changes. This must be taken into account in the case of e.g. colour-coded implants.

neodisher FA can be used with any level of water hardness.

Application and dosage:

neodisher FA is used in washer disinfectors for instruments. The dosage depends on the degree of soiling of the items to be washed and is 2 - 4 ml/l.

Use suitable dosing devices for dosing.

Typical programmes are as follows:

For the automated cleaning and disinfection of surgical instruments and anaesthetic equipment with thermal disinfection:

Cleaning:	2 - 4 ml/l neodisher FA
Neutralisation:	1 - 2 ml/l neodisher Z or neodisher N
Final rinse :	0,3 - 1,0 ml/l neodisher MediKlar

For the removal of blood, serum, foodstuffs and analytical residues in medical and analytical laboratories:

Cleaning:	2 - 4 ml/l neodisher FA
Neutralisation:	1 - 2 ml/l neodisher Z or neodisher N

For the cleaning of baby bottles in milk kitchens:

Cleaning:	2 - 4 ml/l neodisher FA
Neutralisation:	1 - 2 ml/l neodisher Z

For avoiding water stains the use of deionised water in the final rinse is recommended. At the same time this protects anodised aluminium.

The neodisher FA-solution has to be rinsed off completely with water (preferably with deionised water).

Do not mix with other products.

Rinse out dosing system including suction hose with water before changing product. Observe the instrument manufacturer's recommendations for use according to the requirements of DIN EN ISO 17664.

The washer disinfectors' manufacturer's instruction for use are to be observed.

For professional use only.

Technical data:

Density (20 °C): 1.5 g/cm³

pH-range (determined in deionised water, 20 °C): 2 - 4 ml/l: 11.6 - 11.8

Viscosity (concentrate, 20 °C): < 50 mPa s

Conductivity (determined in deionised water, 60 °C): 3 g/l: 3.5 mS/cm

Titration factor: 0.65 (in accordance with neodisher titration method)

Ingredients:

Ingredients according to Regulation (EC) No 648/2004 on detergents:
15 - 30 % phosphates

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With the above information, which is appropriate to our current knowledge we describe our product regarding possible safety necessities, but we do not involve any quality description or promise certain properties.

neodisher[®] FA

Alkaline detergent for the automated reprocessing of thermostable and thermolabile instruments
- liquid concentrate -

CE-marks: CE

neodisher FA conforms to the requirements of the legislation for medical devices according to the European directive 93/42/EC.

Storage information:

Sensitive to frost below - 15 °C. Always store at a temperature range between -15 and 30 °C. Usable for 4 years when stored as recommended. For expiry date refer to the stamp mark on the label behind the hour glass symbol .

Hazard and precautionary statements:

For safety information see EC safety data sheets. These are available at www.drweigert.com under the category "Download".

Dispose only when container is empty and closed. For disposal of product residues, refer to the Material Safety Data Sheet.

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Chemische Fabrik Dr. Weigert GmbH & Co. KG
Chemische Fabrik Dr. Weigert GmbH & Co. KG Telefon: (040) 789 60-0